


>>> DIAQUICK RAPID TESTS



Drugs of Abuse (DOA) in Urine

Single Parameter Tests


Amphetamines (AMP), *Cassette/Dipstick*
Barbiturate (BAR), *Cassette/Dipstick*
Benzodiazepine (BZO), *Cassette/Dipstick*
Buprenorphine (BUP), *Cassette*
Cannabinoids (THC), *Cassette/Dipstick*
Cocaine (COC), *Cassette/Dipstick*
Ethylglucuronide (ETG), *Cassette/Dipstick*
Fentanyl (FYL), *Cassette*
Ketamine (KET), *Cassette*
LSD, *Cassette/Dipstick*
Ecstasy (MDMA), *Cassette/Dipstick*
Methadone (MTD), *Cassette/Dipstick*
Methamphetamine (MET), *Cassette/Dipstick*
Morphine (MOP), *Cassette/Dipstick*
Opiates (OPI), *Cassette*
Oxycodone (OXY), *Cassette*
Phencyclidine (PCP), *Dipstick*
Propoxyphene (PPX), *Cassette*
Spice (K2), *Cassette*
Tramadol (TRA), *Cassette/Dipstick*
Tricyclic Antidepressants (TCA), *Cassette/Dipstick*
Zolpidem (ZOL), *Cassette/Dipstick*

CUT-OFF-VALUES

AMP	1000 ng/mL	MOP	300 ng/mL
BAR	300 ng/mL	MTD	300 ng/mL
BUP	10 ng/mL	OPI	2000 ng/mL
BZO	300 ng/mL	OXY	100 ng/mL
COC	300 ng/mL	PCP	25 ng/mL
ETG	500 ng/mL	PPX	300 ng/mL
FYL	20 ng/mL	Spice	50 ng/mL
KET	1000 ng/mL	TCA	1000 ng/mL
LSD	50 ng/mL	THC	50 ng/mL
MDMA	500 ng/mL	TRA	100 ng/mL
MET	1000 ng/mL	ZOL	50 ng/mL

Multi Parameter Panels

Multi-3 Drug Panel (BZO, COC, MOP)
Multi-3/1 Drug Panel (BUP, MOP, MTD)
Multi-4 Drug Panel (AMP, COC, MOP, THC)
Multi-5 Line Cassette (COC, MET, MOP, MTD, THC)
Multi-5 Drug Panel (COC, MET, BZO, THC, MOP)
Multi-5/3 Drug Panel (AMP, COC, MET, MOP, THC)
Multi-5/4 Drug Panel (AMP, COC, MDMA, MOP, THC)
Multi-5/6 Drug Panel (AMP, BZO, COC, MOP, THC)
Multi-6 Drug Panel (BZO, COC, MET, MOP, MTD, THC)
Multi-6/1 Drug Panel (AMP, BZO, COC, MET, MOP, THC)
Multi-6/4 Drug Panel (AMP, BUP, BZO, MET, MOP, THC)
Multi-6/6 Drug Panel (BUP, COC, MET, MOP, MTD, THC)
Multi-6/7 Drug Panel (BUP, BZO, COC, MOP, MTD, THC)
Multi-7 Drug Panel (AMP, BUP, BZO, COC, MTD, MOP, THC)
Multi-10 Drug Panel
 (AMP, BAR, BZO, COC, MDMA, MET, MOP, MTD, TCA, THC)
Multi-10/1 Drug Panel
 (AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, THC)
Multi-10/2 Drug Panel
 (AMP, BAR, BZO, COC, MDMA, MOP, MTD, OPI, TCA, THC)
Multi-10/3 Drug Panel
 (AMP, BZO, COC, MDMA, MOP, MTD, OPI, PCP, TCA, THC)
Multi-10/4 Drug Panel
 (AMP, BAR, BUP, BZO, COC, MDMA, MET, MTD, OPI, THC)
Multi-10/5 Drug Panel
 (AMP, BAR, BUP, BZO, COC, MET, MOP, MTD, TCA, THC)
Multi-10/6 Drug Panel
 (AMP, BAR, BZO, COC, MET, MOP, MTD, PCP, TCA, THC)
Multi-10/7 Drug Panel
 (AMP, BAR, BZO, COC, MET, MTD, OPI, PCP, TCA, THC)
DOA Multi-6 Cup (AMP, BZO, COC, MET, MOP, THC)
Multi-6 drug testing panel combined with cup and urine temperature control



Urine Adulteration

Drug Adulteration Strip, *colorimetric test*

DOA Saliva **CE**

Alcohol Saliva Strip

DOA Saliva Multi-6 Test

(AMP, COC, MET, OPI, PCP, THC)

DOA Saliva Multi-6/1 Test

(AMP, COC, MET, OPI, MTD, THC)

both also available as High Sensitivity Touch Test

DOA Saliva Twist Test

(AMP, COC, MET, MTD, OPI, THC)

Pregnancy / Fertility **CE**

hCG Dipstick *(urine)*

hCG Combo Dipstick *(urine/serum)*

hCG Cassette *(urine)*

hCG Combo Cassette *(urine/serum)*

hCG Midstream Test *(urine)*

LH Dipstick *(urine)*

LH Cassette *(urine)*

Tests for Home Use **CE**

Pregnancy Test - Dipstick *(urine)*

DiaFem Pregnancy Test - Midstream *(urine)*

Fertility Test - Dipstick *(urine)*

Fertility Test - Cassette *(urine)*

DiaFem Fertility Test - Midstream *(urine)*

Tumour Markers **CE**

AFP Cassette *(serum/plasma/whole blood)*

CEA Cassette *(serum/plasma/whole blood)*

FOB Cassette *(stool samples)*

FOB-TRF Combi Cassette *(stool samples)*

PSA Cassette *(serum/plasma)*

Cardiac Markers **CE**

Cardiac Combo Cassette *(serum/plasma/whole blood)*

Troponin I Cassette *(serum/plasma/whole blood)*

Others **CE**


CRP Cassette *(serum/plasma/whole blood)*

D-Dimer Cassette *(plasma/whole blood)*

Microalbumin Dipstick *(urine)*



DIAQUICK RAPID TESTS



Aids, Hepatitis

HAV IgM Cassette (serum/plasma/whole blood)

HBsAb Cassette (serum/plasma)

HBsAb Dipstick (serum/plasma)

HBsAg Cassette (serum/plasma) **CE**

HBsAg Cassette (serum/plasma/whole blood)

HBsAg Dipstick (serum/plasma)

HCV Cassette (serum/plasma) **CE**

HCV Cassette (serum/plasma/whole blood)

HIV Cassette (serum/plasma)

HIV Cassette (serum/plasma/whole blood)

HIV Triline Cassette (serum/plasma/whole blood) **CE**

ToRCH **CE**

Toxoplasma IgM Cassette (serum/plasma)

Toxoplasma IgG/IgM Cassette (serum/plasma)*

Rubella IgM Cassette (serum/plasma)

CMV IgM Cassette (serum/plasma)

HSV 1/2 IgM Cassette (serum/plasma)

ToRCH IgM Combo Cassette (serum/plasma)

* no **CE**

Other Infectious Diseases **CE**

Chlamydia Cassette (endocervical/urethral/urine samples)

Dengue IgG/IgM Ab Cassette (serum/plasma/whole blood)

H.pylori Cassette (serum/plasma/whole blood)

H.pylori Stool Cassette (stool samples)

Influenza Ag Dipstick (nasal/throat swab samples)

Malaria P.f. Cassette (whole blood)

Malaria P.f./Pan Cassette (whole blood)

Malaria P.f./P.v. Cassette (whole blood)

Mononucleosis Cassette (serum/plasma/whole blood)

Adenovirus Cassette (stool samples)

Rotavirus Cassette (stool samples)

Rota/Adeno Combi Cassette (stool samples)

S.typhi IgG/IgM Ab Cassette (serum/plasma/whole blood)

S.typhi Ag Cassette (stool/serum/plasma/whole blood)

S.typhi/paratyphi Ag Cassette (stool/serum/plasma/whole blood)

Strep.A Cassette (throat swab samples)

Strep.A Dipstick (throat swab samples)

Syphilis Cassette (serum/plasma/whole blood)

Syphilis Dipstick (serum/plasma/whole blood)